


Tips to make it stick better:

- clean and dry skin
- the less hair the better

If it doesn't stick perfectly the first attempt:


- test another patch on another area of the body
- 96.9% of 'sticking' problems are due to process, not the adhesive

1



Peel off one side.

2



Align middle 'non-stick' section over the sensor*.


3



Gently smooth down the adhesive to the skin and lay over the sensor*.


4

5




Don't touch exposed adhesive.
Touching the adhesive reduces stick-ability.

6



Peel back remaining side and smooth over the skin.

7



Gently rubbing the patch helps activate the adhesive.

8

Have a splendid day...
and **STICK** to your dreams!

**Supports Libre and Dexcom G6 sensors.*

Want to give your NJaP Patch the best chance of success?

Keep reading!

What you do before and after applying NJaP matters....trust us...we know.

1 Your skin must be clean and dry and free of oily things like moisturisers and the less hairy the area the better sticking ability your patch will have.

2 Avoid placing the patch near creases or fold in the skin - for example if you wear your CGM high on the back of your arm and the patch edge is close to the underarm area, the movement in this area will likely cause the patch to come loose in that area.

3 The less active you are after applying the patch the better chance the NJaP Patch has of bonding with your skin. Ideally, no exercise or sweating for 2 hours and even longer if possible. This bonding period is crucial to the next 10 days of patch adhesive ability. If you sweat within 2 hours after patch application, your patch is unlikely to last more than a few days. Definitely do not shower within a few hours.

4 Clothes rubbing on the patch within a few hours of application will also cause the patches to last shorter periods. Ideally no interaction with the patch within a few hours to ensure best sticking ability.

5 Avoid stretching the patch on application and try not to touch the adhesive area during application, this reduces sticking ability.

6 NJaP Patches come with a non-stick backing section. We do this to give you the choice of whether you want the patch to stick to your device. Not sticking to the device allows you to change your patch as needed without the need to remove your device. If you prefer the patch to stick to your device, by all means remove the middle section and place your patch on top of the device.

7 If after following these tips you find your patch is coming loose, we recommend that you try the patch on another section of skin that is clean and dry and free of hair. NJaP Patches have a highly effective hypoallergenic latex free adhesive substance that is consistent across all patches. We have tested countless patches that have been claimed to not stick and 100% of these have worked on other people, for this reason we know that it's always the skin and the application process that determines the patch's success.

8 If you think the patches aren't working and staying stuck on, we recommend trying a patch on a different area of your body with no sensor. This will help confirm if it is the adhesive or your skin or the area your CGM is on your body.

9 We highly recommend using a product like Skin Glu which helps to add a barrier on your skin and also gives your skin more stickiness. Best to apply Skin Glu after you have used Skin Glee, which is an alcohol wipe designed to ensure the cleanest skin for device application.

10 Finally, there will be a small percentage whose skin for any number of reasons just isn't friendly to having things stick to it e.g. maybe naturally oily skin - we estimate this is 2 to 3% of the population.

**NOT
JUST A
PATCH™**